



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Volumetric pump ALARIS CAREFUSION GW (Reconditioned)

Used,
Technical condition: very good,
Visual condition: very good,
Manufacture U.K., (2012),
Power supply: 230 V,
Emergency power supply: through 1 built-in rechargeable NiMH 7.2 V 1.25 Ah battery
Frequency: 50/60 Hz,
Power: 10 VA,
Alaris GW is a small, lightweight volumetric infusion pump for accurate and reliable infusions,
Can be used for intravenous, intra-arterial, subcutaneous and epidural administration of fluids and drugs,
Designed for use in general and intensive care,
In-line air detection,
Foiled control panel,
German language keyboard on control panel,
Built-in battery ensures uninterrupted pump operation in case of power loss or patient transport,
Housing without protruding knobs and buttons on the front panel for easy maintenance and disinfection,
Can be mounted on a tripod,
Safety features:
BackOff function to reduce post occlusion bolus,
Manual shut-off,
Accumulated air detection in line,
Automatic recognition of intravenous infusion sets,
Alarm states:
System error,
Air in line,
Open door,
Invalid infusion line,
Door open,
Battery discharged,
Upper part occlusion,
Occlusion in lower part,
Equipped with:
RS - 232 socket for nurse call button,
Socket for flow sensor,
Holder with clamping screw for mounting on a rail, e.g. stand or bed,
LED screen,
Specifications:
Infusion rate range:
Standard mode: 1 - 999 mL / h,
Micro mode: 1 - 99.9 mL / h,
Volume to be administered (VTBI):
Standard mode: 1 - 9999 mL / h,
Micro mode: 0.1 - 999 mL / h,
Volume administered (VI):
Standard and Micro modes: 0 - 9999 mL,
KVO infusion rate: max 5 mL / h or equal to infusion rate,
Occlusion alarm pressure (at 125 mL / h rate):
Low: 250 mmHg,
Normal: 350 mmHg,
High: 500 mmHg,
Bolus administration during infusion:
Bolus delivery rate: 1- 999 mL / h,
Bolus volume: 0 - 99 mL,
Included:
Ivac 180 flow sensor,
Operating instructions in languages: PL, EN, DE, CS, BG, ES, FR, IT, NL, PT, ZH (PDF),
Power cord,
Dimensions: 13.5 x 10.5 x 13 cm,
Weight: 2 kg,
Cat 02,
Has a current Technical Passport issued,
Warranty:
6 months for the domestic market - Poland,
1 month for the international market,
Financing options (Poland only): Installments, Leasing, Loan,
If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.
WE SPEAK ENGLISH
HABLAMOS ESPANOL
WIR SPRECHEN DEUTSCH
ON PARLE FRANCAIS
WIJ SPREKEN NEDERLANDS
МЫ ГОВОРИМ ПО РУССКИ
МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ
WE SHIP WORLDWIDE

