



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Infusion pump docking station FRESENIUS VIAL ORCHESTRA BASE A (Reconditioned)

Used,
Technical condition: good,
Visual condition: good,
French production,
Power supply: 230 V,
Emergency power supply: through 1 new built-in 6 V 1.1/1.3 Ah battery
Frequency: 50/60 Hz,
Power: 125 VA,
The Base A docking station is a data hub and power supply for up to 6 DPS modules and/or DPS Visio modules and/or MVP modules,
Allows faster therapeutic control,
Communication interface compatible with the PDMS main software,
Built-in RS-232 interface,
Indicator lights for power supply and infusion pump operation,
No pumps included - can be retrofitted with FRESENIUS ORCHESTRA infusion pumps (syringe or drip) at additional cost,
Included:
Operating instructions in German (PDF),
Power cord,
Dimensions: 33.5 x 26 x 8 cm,
Weight: 2 kg,
Cat 03,
Has a current Technical Passport issued,
Warranty:
6 months for the domestic market - Poland,
1 month for the international market,
Financing options (Poland only): Installments, Leasing, Loan,
If you have any questions, please do not hesitate to contact us!
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.
WE SPEAK ENGLISH
HABLAMOS ESPANOL
WIR SPRECHEN DEUTSCH
ON PARLE FRANCAIS
WIJ SPREKEN NEDERLANDS
МЫ ГОВОРИМ ПО РУССКИ
МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ
WE SHIP WORLDWIDE

