



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Infusion pump docking station FRESENIUS VIAL ORCHESTRA BASE A (Reconditioned)**

Used,  
Technical condition: very good,  
Visual condition: very good,  
French production,  
Power supply: 230 V,  
Emergency power supply: 1 new built-in 6 V 1.1/1.3 Ah battery  
Frequency: 50/60 Hz,  
Power: 125 VA,  
The Base A docking station is a data hub and power supply for up to 6 DPS modules and/or DPS Visio modules and/or MVP modules,  
Allows faster therapeutic control,  
Communication interface compatible with the PDMS main software,  
Built-in RS-232 interface,  
Indicator lights for power supply and infusion pump operation,  
No pumps included - can be retrofitted with FRESENIUS ORCHESTRA infusion pumps (syringe or drip) at additional cost,  
Included:  
Operating instructions in German (PDF),  
Power cord,  
Dimensions: 33.5 x 26 x 8 cm,  
Weight: 2 kg,  
Cat 04,  
Has a current Technical Passport issued,  
Warranty:  
6 months for the domestic market - Poland,  
1 month for the international market,  
Financing options (Poland only): Installments, Leasing, Loan,  
If you have any questions, please do not hesitate to contact us!  
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.  
WE SPEAK ENGLISH  
HABLAMOS ESPANOL  
WIR SPRECHEN DEUTSCH  
ON PARLE FRANCAIS  
WIJ SPREKEN NEDERLANDS  
МЫ ГОВОРИМ ПО РУССКИ  
МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ  
WE SHIP WORLDWIDE

