



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



OLYMPUS FORCE MAX NM - 401L - 0625 injection needle (Ex-Display)

New, demo,
Japanese production,
Force Max Olympus is a sterile injection needle for a wide range of endoscopic applications,
Optimal for upper gastrointestinal tissues,
Compatible with Olympus endoscopes: EF, GIF, OGF, CF (except I and L length), OSF,
Extremely rigid bend-resistant tube allows for deeper penetration,
A locking mechanism with a clearly audible click indicates when the needle blade is fully retracted into the sheath,
Equipped with:

Handle with slider for inserting and withdrawing the needle,
fluid supply tube,
Needle for puncturing the mucosa for vascular or submucosal injection,
Port for syringe connection and drug administration,

Specifications:

Needle length: 6 mm,
Needle diameter: 0.5 mm (25G),
Needle blade bevel angle: 14°,
Working length of tool: 1650 mm,
Minimum working channel diameter: 2 mm,
Maximum diameter of the part introduced into the endoscope: 1.9 mm,

Included:

Operating instructions in EN, DE languages (PDF),

Handle dimensions: Ø 1 x 16 cm,

Weight of set: 0.5 kg,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

