



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



OLYMPUS NM - 201L - 0625 injection needle (Ex-Display)

New, demo,
Japanese production,
Olympus NM - 201L is a sterile injection needle with a wide range of applications in endoscopy,
Optimal for upper gastrointestinal tissues,
Equipped with:

Handle with slider for insertion and withdrawal of the needle,
Delivery tube,
Needle for puncturing the mucosa for vascular or submucosal injection,
Port for syringe connection and drug administration,

Specifications:

Needle length: 6 mm,
Needle diameter: 0.5 mm (25G),
Needle blade bevel angle: 14°,
Working length of tool: 1650 mm,
Minimum working channel diameter: 2 mm,
Maximum diameter of the part introduced into the endoscope: 1.9 mm,
Handle dimensions: 13.5 x 2 x 1 cm,
Weight: 0,5 kg,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

