



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



OLYMPUS ENDOJAW FB - 210U biopsy forceps (Ex-Display)

New, demo,
Japanese production,
Endojaw Olympus are fenestrated biopsy forceps of the "alligator" type adapted to tissue collection,
Optimal for lower gastrointestinal tissues,
They have an intuitive to use handle with a slider to open and close the forceps,
The sheath design ensures smooth insertion of the forceps and the two-stage blade ensures the precision of the procedure,
The pendulum design of the forceps jaws allows for easier positioning and better control of the forceps,
Low risk of canal damage,

Specifications:

Tool working length: 2300 mm,
Minimum working channel diameter: 2.8 mm,
Maximum diameter of the part entering the endoscope: 2.45 mm,
Handle dimensions: 16 x 3.5 x 3.5 cm,
Weight: 0.5 kg,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

