



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Docking station for 6 BECTON pumps - DICKINSON BASE A (Reconditioned)

Used,
Technical condition: very good,
Visual condition: very good,
French production,
Power supply: 230 V,
Emergency power supply: 1 new built-in 6 V 1.3 Ah battery
Frequency: 50/60 Hz,
Power: 55 VA,
The Base A docking station is a data hub and power supply for up to 6 DPS modules and/or DPS Visio modules and/or MVP modules,
Allows faster therapeutic control,
Communication interface compatible with the PDMS main software,
Built-in RS-232 interface,
Indicator lights for power supply and infusion pump operation,
One single power cable for all connected modules,
Equipped with stand with convenient carrying handle,
No pumps included - can be retrofitted with FRESENIUS ORCHESTRA infusion pumps (syringe or drip) at additional cost,
Included:
Operating instructions in German (PDF),
Power cord,
Dimensions:
General: 33.5 x 29 x 8 cm,
Station height with tripod: 80 cm,
Weight of station:
Without pumps: 3.5 kg,
With pumps: 17.5 kg,
It has a current Technical Passport issued,
Warranty:
6 months for the domestic market - Poland,
1 month for the international market,
Financing options (Poland only): Installments, Leasing, Loan,
If you have any questions, please do not hesitate to contact us!
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.
WE SPEAK ENGLISH
HABLAMOS ESPANOL
WIR SPRECHEN DEUTSCH
ON PARLE FRANCAIS
WIJ SPREKEN NEDERLANDS
МЫ ГОВОРИМ ПО РУССКИ
МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ
WE SHIP WORLDWIDE

