



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### **Medical chair (Ex-Display)**

New, demo,  
Technical condition: perfect,  
Visual condition: perfect,  
Medical chair is a mobile, elegant furniture with a stable construction designed for doctors surgeries, clinics and hospitals,  
Upholstery with medical certificate in green (possibility to exchange the upholstery for a new one from our colour chart at extra cost),  
Equipped with comfortable armrests,  
Seat height adjustable from 41 to 54 cm with a gas spring,  
Backrest height adjustable within 108 - 114 cm,  
Backrest can be tilted in the range 0 - 35°,  
Mobile on a five-star base with double castors,  
Dimensions:  
Seat: 46 x 42.5 x 5 cm,  
General: Ø 60 x 108/114 cm,  
Weight: 12 kg,  
Comes with a current Technical Passport issued,  
Warranty:

6 months for the domestic market - Poland,

1 month for the international market,

Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

