



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Sinusoscope ZHEJIANG TIANSONG MEDICAL INSTRUMENT B0310 (New)

New,
P.C.R. production,
Tiansong Sinusoscope is a rigid nasal endoscope for sinus examination and surgery designed for ENT offices and clinics,
External diameter: Ø 2.7 mm,
Working length: 175 mm,
Observation angle: 0°,
Can be sterilized in autoclave at 2 bar - 5 min - 134° C,
Optics compatible with Storz and Wolf standards,
Access to Wolf connector by unscrewing the Storz tip,
Dimensions: Ø 3 x 23.5 x 5.5 cm,
Weight: 0.5 kg,
Holds a valid technical passport,
Warranty: 12 months,
Financing options (Poland only): Instalments, Leasing, Loan,
If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

