



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Unit combi ZIMMER SONOMAT II, SINUS, VACO (Reconditioned)

Used,
Technical condition: very good,
Visual condition: very good,
German production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Kombi Zimmer is a set of devices intended for physical therapy procedures,
The combo set consists of:
ZIMMER ultrasound therapy device Model SONOMAT II:
Power: 60 W,
The therapy is based on the action of ultrasound waves on deeper structures in order to achieve therapeutic effects, e.g. analgesic,
Adjustable intensity and frequency,
Electrostimulator ZIMMER Model SINUS:
Power: 45 VA,
Electrostimulation uses low and medium frequency current to reduce pain, reduce swelling, relax muscles and improve conductivity,
ZIMMER Vacuum Therapy Device Model VACO:
Power: 80 VA,
Has a built-in manometer,
Keypad on control panel in German,
All placed on a mobile trolley on 4 double wheels,
Included:
ERBE vacuum therapy electrodes with Ø 6 cm - 4 pcs,
Tube for connecting the vacuum electrodes,
Cable for connecting the electrostimulator,
ZIMMER ultrasound head with a frequency of 795 kHz and an area of 5 cm²,
Cable for connecting two devices,
Overall dimensions: 61 x 47 x 89 cm,
Weight: 41 kg,
It has a current Technical Passport issued,
Warranty:
6 months for the domestic market - Poland,
1 month for the international market,
Financing options (Poland only): Installments, Leasing, Loan,
If you have any questions, please do not hesitate to contact us!
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.
WE SPEAK ENGLISH
HABLAMOS ESPANOL
WIR SPRECHEN DEUTSCH
ON PARLE FRANCAIS
WIJ SPREKEN NEDERLANDS
МЫ ГОВОРИМ ПО РУССКИ
МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ
WE SHIP WORLDWIDE

