



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Patient lift with hygienic cut-out (bathtub) ARJO MECANAIDS AMBULIFT B (Reconditioned)

Used,
Technical condition: good,
Visual condition: very good,
Dutch manufacture,
The Mecanaids Ambulift is a robust bath lift designed for safe and comfortable lifting and transfer of patients during bathing,
With its intuitive operation and sturdy design, it is suitable for use in the home, hospital wards, rehabilitation centres and nursing homes,
The seat has a hygienic cut-out, which makes it easier to carry out nursing activities,
smooth manual seat height adjustment from 17 to 104 cm with a crank,
The seat can be removed,
Equipped with 2 patient handles,
Mobile, on a mobile base on 4 castors Ø 8 cm, 2 of which with brakes,
Dimensions:
Seat: 43.5 x 35 cm,
Backrest height: 38 cm,
Overall: 100.5 x 78 x 153 cm,
Maximum load: 160 kg,
Weight: 48 kg,
Comes with a current Technical Passport issued,
Warranty:
6 months for domestic market - Poland,
1 month for the international market,
Financing options (Poland only): Installments, Leasing, Loan,
If you have any questions, please do not hesitate to contact us!
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.
WE SPEAK ENGLISH
HABLAMOS ESPANOL
WIR SPRECHEN DEUTSCH
ON PARLE FRANCAIS
WIJ SPREKEN NEDERLANDS
МЫ ГОВОРИМ ПО РУССКИ
МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ
WE SHIP WORLDWIDE

