



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MEDRAD MARK V PRO VIS PPD angiography syringe (Reconditioned)

Used,
Technical condition: very good,
Visual condition: very good,
U.S.A. production,
Frequency: 50/60 Hz,
The Medrad Mark V Pro Vis is a syringe for contrast administration for angiographic studies,
It is designed to be easy to use, from the legibility of the display and controls to the ergonomics and functionality of the injector and control panel,
The programmed microprocessor helps protect against excessive volume, flow and pressure values,
Constantly updated syringe volume,
Syringe seated on a two-piece arm allowing its height to be adjusted from 131 - 159 cm,
Possibility of rotating the control panel by 180°,
Equipped with a mechanical stop that is automatically set and locked to physically limit the injection to the selected volume,
The set position is verified by LEDs located on the injector head,
Up to four flow rates can be programmed,
Easy integration into any imaging system,
Control panel divided into sections for efficient syringe programming,
Large, secure handles and manoeuvrable castors make the unit easy to manoeuvre,
Clear and easy to read control panel is easy to keep clean due to foiling,
Keypad on control panel in language: English,
Specifications:
Charging speed: 5 - 10 mL/s,
Flow speed:
For 60 mL syringe:
0.3 - 10 mL/s or mL/min or mL/h, Step 0.1 mL,
10 - 20 mL/s, step 1 mL,
10 - 59 mL/min or mL/h, step 1 mL,
For 150 and 200 mL syringes:
0.3 - 10 mL/s or mL/min or mL/h, step 0.1 mL,
10 - 50 mL/s, step 1 mL,
10 - 59 mL/min or mL/h, step 1 mL,
Volume: from 0.5 mL for 60 mL syringe or from 1 mL for 150 and 200 mL syringes up to the maximum volume of the selected syringe, step 0.1 mL,
Maximum programmable pressure limit: 7584 kPa for mL/s unit,
Scanning delay: 0 - 99.9 s, step 0.1 s,
Injection time: 0 - 650 s,
Flow rate rise/fall time: 0 - 9.9 s, step 0.1 s,
Syringe heating: 37° C,
Syringes supported: 60, 150, 200 mL,
Mobile, on a base with Ø 12 cm castors, all with brake,
Included:
MEDRAD pressure jackets for 60 mL syringe - 2 pcs,
Operating instructions in Spanish (PDF),
Overall dimensions: 74 x 71 x 131/159 cm,
Weight: 97.5 kg,
Has a current Technical Passport issued,
Warranty:
6 months for domestic market - Poland,
1 month for the international market,
Financing options (Poland only): Installments, Leasing, Loan,
If you have any questions, please do not hesitate to contact us!
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.
WE SPEAK ENGLISH
HABLAMOS ESPANOL
WIR SPRECHEN DEUTSCH
ON PARLE FRANCAIS
WIJ SPREKEN NEDERLANDS
МЫ ГОВОРИМ ПО РУССКИ
МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ
WE SHIP WORLDWIDE

