



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Video - gastroscope OLYMPUS EVIS EXERA II GIF - Q165 (Reconditioned)**

Used,  
Technical condition: very good,  
Visual condition: very good,  
Japanese production,  
Olympus Evis Exera II is a video gastroscope that allows for high quality images,  
The outer diameter of the endoscopic probe tip and the outer diameter of the endoscopic probe tube provide excellent insertion capabilities,  
Scope ID function stores information in the built-in memory chip and displays it on the monitor,  
Wide field of view allows accurate observation of the selected area,  
The 2.8 mm diameter endoscope channel allows the use of a wide range of therapy accessories,  
Convenient access to buttons and user programmable switches for easy operation,  
Endoscope tip bendable in 4 directions allows complete and comprehensive examination of the upper gastrointestinal tract,  
Technical parameters:  
Optical system:  
Field of view: 140°,  
Viewing direction: forward view,  
Depth of focus: 3 - 100 mm,  
Endoscopic probe:  
Outer diameter of endoscopic probe tip: 9.2 mm,  
Outer diameter of endoscopic probe tube: 9.2 mm,  
Working length: 1030 mm,  
Biopsy channel:  
Inside diameter of the channel: 2.8 mm,  
Minimum viewing distance: 3.0 mm from the tip of the endoscopic probe,  
Tip bendable in range:  
210° up,  
90° down,  
100° right,  
100° left,  
Total length: 1345 mm,  
Compatible with system: CV -160 / 145 / 140,  
Included:  
OLYMPUS MH - 443 suction valve,  
Water and air valve OLYMPUS MH - 438,  
Lid,  
Biopsy valve OLYMPUS VB - 301,  
Biopsy valve,  
OLYMPUS carrying case,  
Case dimensions: 60 x 12 x 43 cm,  
Weight: 5 kg,  
Cat. 08,  
It has a current Technical Passport issued,  
Warranty:  
6 months for the domestic market - Poland,  
1 month for the international market,  
Financing options (Poland only): Installments, Leasing, Loan,  
If you have any questions, please do not hesitate to contact us!  
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.  
WE SPEAK ENGLISH  
HABLAMOS ESPANOL  
WIR SPRECHEN DEUTSCH  
ON PARLE FRANCAIS  
WIJ SPREKEN NEDERLANDS  
МЫ ГОВОРИМ ПО РУССКИ  
МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ  
WE SHIP WORLDWIDE

