



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



BAYER MEDRAD STELLANT CT D Syringe for contrast injection (Reconditioned)

Used,

Technical condition: very good,

Visual condition: very good,

U.S.A. production,

Power supply: 230 V,

Frequency: 50/60 Hz,

Power: 300 VA,

Medrad Stellant is a precision, dual-headed syringe for contrast administration offering real-time injection graphics,

The syringe consists of:

Syringe head,

11" monitor,

Base unit,

Product features:

Ability to inject simultaneously from both cartridges (contrast and saline) with programming of the percentage of solution injected,

Provides saline flushing capability and extensive protocol options for the most advanced clinical applications,

Automatic piston movement and syringe filling,

Allows variable proportions of syringe plunger movement with contrast and saline simultaneously,

Integrated saline test injection,

Capable of test and multiphase injection up to 6 phases / protocol,

Heating of cartridges (contrast and saline),

Designed for complex CTA and cardiac CT protocols,

Creates a compact contrast bolus for better visualization,

With proper ratio left and right ventricle can be illuminated evenly improving image quality,

It has two digital displays,

Keypad on control panel in English,

Can store up to 32 injection protocols,

Allows graphical pressure monitoring during injection,

Specifications:

Flow rate: 0.1 - 10 mL, step 0.1 mL,

Volume: from 1 mL to maximum syringe volume, step 1 mL,

Programmable pressure limit: 345, 689, 1034, 1379, 1551, 1724, 2068, 2241 kPa,

Scan delay: 0 - 300 s, step 1 s,

Pause: 0 - 900 s, step 1 s,

Injection time: 0 - 999 s,

Syringe heating: 35° C ± 5° C,

Syringe support: 200 mL,

Mobile, on a trolley with 5-arm MEDRAD Model SHP 800:

Has double castors with Ø 8 cm, 2 of them with brake,

Equipped with a hanger with two fluid hooks,

Included:

2 communication links,

Instruction manual in PL, EN (PDF),

Power cord,

Dimensions:

Syringes on trolley: 78 x 65 x 142 cm,

Monitor: 32 x 22 x 34 cm,

Base: 28 x 22 x 27.5 cm,

Weight of set: 53 kg,

Comes with a current Technical Passport issued,

Warranty:

6 months for the domestic market - Poland,

1 month for the international market,

Financing options (Poland only): instalments, leasing, loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

