



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### **Certified rapid cassette test for SARS-COV-2 antigen (saliva) REALY TECH 1 pcs pack (New)**

New,  
P.C.R. production,  
Realy Tech is an in vitro diagnostic cassette test for the qualitative detection of SARS - Cov - 2 antigens in a saliva sample using a rapid immunochromatographic method,  
Rapid and easy collection of saliva samples reduces the time needed to perform the test,  
Identification is based on the detection of monoclonal antibodies that are specific for coronavirus antigen,  
Reading of the result possible after 10 - 20 minutes,  
Easy and self-sampling reduces the risk of virus transmission during the test,  
Avoids problematic throat or nasopharyngeal sampling, making the test also suitable for children,  
Comes with CE declaration of conformity and can be marketed in the Republic of Poland,  
Specification:  
Sensitivity: 92.9 %,  
Specificity: 99.58 %,  
Included:  
Test cassette,  
Collector for saliva collection,  
Extraction tube with stopper and dropper,  
Ampoule with buffer,  
Leaflet in EN, DE,  
If you have any questions, please do not hesitate to contact us!  
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.  
WE SPEAK ENGLISH  
HABLAMOS ESPANOL  
WIR SPRECHEN DEUTSCH  
ON PARLE FRANCAIS  
WIJ SPREKEN NEDERLANDS  
МЫ ГОВОРИМ ПО РУССКИ  
МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ  
WE SHIP WORLDWIDE

