



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Light source RICHARD WOLF 4015 (Reconditioned)

Used,
Technical condition: very good,
Visual condition: very good,
German production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 220 VA,
Wolf 4015 is a halogen light source for diagnostic and surgical procedures,
Equipped with 2 bulbs 150 W / 24 V,
Compatible with light guides of Ø 8 - 11 mm,
Infinitely adjustable light intensity,
A knob on the front of the unit allows you to change between the main and emergency bulbs and to release the lock to access them,
Comes with carrying handle,
Power cord included,
Dimensions: 26.5 x 19.5 x 14.5 cm,
Weight: 5.5 kg,
Cat 03,
Has a current Technical Passport issued,
Warranty:
6 months for the domestic market - Poland,
1 month for the international market,
Financing options (Poland only): Installments, Leasing, Loan,
If you have any questions, please do not hesitate to contact us!
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.
WE SPEAK ENGLISH
HABLAMOS ESPANOL
WIR SPRECHEN DEUTSCH
ON PARLE FRANCAIS
WIJ SPREKEN NEDERLANDS
МЫ ГОВОРИМ ПО РУССКИ
МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ
WE SHIP WORLDWIDE

