



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### **Light source RICHARD WOLF 4015 (Reconditioned)**

Used,  
Technical condition: very good,  
Visual condition: very good,  
German production,  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Power: 220 VA,  
Wolf 4015 is a halogen light source for diagnostic and surgical procedures,  
Equipped with 2 bulbs 150 W / 24 V,  
Compatible with light guides of Ø 8 - 11 mm,  
Infinitely adjustable light intensity,  
A knob on the front of the unit allows you to change between the main and emergency bulbs and to release the lock to access them,  
Comes with carrying handle,  
Power cord included,  
Dimensions: 26.5 x 19.5 x 14.5 cm,  
Weight: 5.5 kg,  
Cat 03,  
Has a current Technical Passport issued,  
Warranty:  
6 months for the domestic market - Poland,  
1 month for the international market,  
Financing options (Poland only): Installments, Leasing, Loan,  
If you have any questions, please do not hesitate to contact us!  
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.  
WE SPEAK ENGLISH  
HABLAMOS ESPANOL  
WIR SPRECHEN DEUTSCH  
ON PARLE FRANCAIS  
WIJ SPREKEN NEDERLANDS  
МЫ ГОВОРИМ ПО РУССКИ  
МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ  
WE SHIP WORLDWIDE

