



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **CU MEDICAL SYSTEM MEPAD CU - SP1 defibrillator (New)**

- New,
- Manufactured in Korea,
- Power supply: 230 V,
- Emergency power supply: 1 replaceable 12 V 4.2 Ah battery,
- Frequency: 50/60 Hz,
- **The Me Pad CU-SP1** is a semi-automatic defibrillator specifically designed for use by people who are not doctors or paramedics,
- Its safe use, compact size and light weight make it suitable for use in public areas,
- The device guides the assisting person through the CPR process with clear and understandable voice commands and indicators - LED and graphic,
- Equipped with two electrode sizes, making it suitable for both adults and children,
- The possibility to activate the paediatric mode without having to change the electrodes speeds up the assistance,
- The defibrillator performs internal monitoring of electrode quality,
- Features a built-in self-tests,
- Function of automatic adjustment of alarm volume to the environment,
- **Voice commands in PL, EN, DE,**
- The device saves all relevant data from the last 5 resuscitations,
- Battery provides operation for 5 years or 200 shocks,
- Can be transported safely thanks to a dedicated bag (**bag not included, can be retrofitted at additional cost**),
- **Specifications:**
  - Output energy at 50  $\Omega$ :
  - Adult patients: 150 J,
  - Paediatric patients: 50 J,
  - Shock charging time:  $\leq 8$  s,
- **The set includes:**
  - Adult electrode set,
  - Pediatric electrode set,
  - User manual in German,
- Dimensions: 26 x 25.5 x 7 cm,
- Weight: 2.4 kg,
- Valid technical passport issued,
- Warranty 12 months,
- Financing available (Poland only): Instalments, Leasing, Loan,

