



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



AESCULAP video processor PV150 with light source OP 927 and camera head PV151 (set) (Reconditioned)

Used,
Technical condition: very good,
Visual condition: very good,
German production,
The set includes:
AESCULAP Video Processor Model DAVID DIGITAL PV150:
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 25 VA,
Intuitive operation,
Compatible with rigid endoscopes,
White balance control,
Foiled control panel,
Has the following ports: S-VIDEO OUT, VIDEO OUT, SDI OUT, REMOTE,
Foiled control panel,
Keypad on control panel in English,
Dimensions: 30.5 x 30.5 x 8 cm,
Weight: 4 kg,
Light source AESCULAP Model 300 OP 927:
Power supply: 230 V,
Frequency: 50/60 Hz,
Compatible with rigid endoscopes,
Foiled control panel,
Main and emergency lighting: 24 V / 250 W halogen bulb - 2 pcs,
8-step regulation of light intensity,
Colour temperature: approx. 3400 K,
Foiled control panel,
Dimensions: 30.5 x 30.5 x 13 cm,
Weight: 8.5 kg,
AESCULAP Camera Head Model PV151:
Equipped with two buttons for easy control,
Power cord - 2 pcs,
Weight of set: 13 kg,
It has a current Technical Passport issued,
Warranty:
6 months for domestic market - Poland,
1 month for the international market,
Financing options (Poland only): instalments, leasing, loan,
If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

