



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Medical laser QUANTA SYSTEM Q1 QUANTA C* (Reconditioned)

Used,

Technical condition: very good,

Visual condition: very good,

Made in Italy, (2010),

Power supply: 230 V,

Frequency: 50/60 Hz,

Power: 300 VA,

Quanta C* is a combined 980 nm IR diode laser (Quanta C) and 532 nm green laser (Quanta D),

The foot attachment and the carrying case bear signs of use,

Device equipped with:

Built-in 6" touch screen,

Safety button "STOP",

Calibration port,

Key lock to lock the unit,

Remote door lock,

Choice of menu languages: EN, IT,

2 ports for fibre optics connection,

RS - 232 port,

Specification:

Laser type: diode laser,

Modes of operation: CW, pulse, pulse series,

Beam delivery method: optical fibre,

Wavelength (nm):

532 nm,

980 nm,

Pulse length: 3 - 2500 ms,

Repetition rate: 1 - 200 Hz,

Treatment spot size: 0.6 / 0.8 / 1.2 / 1.8 / 2.4 mm,

Number of pulses in series: 1 - 50,

Timer: 1 s - 10 min,

Included:

Head with optical fibre - 2 pcs,

Distance pieces - 6 pcs,

Remote door lock,

Foot attachment,

Key for the lock - 2 pcs,

Instruction manual in Polish,

Suitcase,

Power cord,

Dimensions: 39,5 x 39 x 25 cm,

Weight: 16.5 kg,

Comes with a current Technical Passport issued,

Warranty:

6 months for the domestic market - Poland,

1 month for the international market,

Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

