



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



**ERBE VIO 200D surgical diathermy with APC2 module on trolley (Reconditioned)**

Used,  
Technical condition: very good,  
Visual condition: very good,  
German production,  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Power: 500 W / 920 VA,  
Erbe Vio 200D is a surgical diathermy applicable in gastroenterology, endoscopy and pulmonology,  
System features and benefits:  
Modular design,  
Colourful, easy to read display,  
Selectable menu in the following languages: PL, EN, DE, RU, CS,  
Can be individually configured according to operator needs  
Voltage regulation for gentle, reproducible cutting and coagulation,  
Electric arc control for high energy cutting or coagulation and for underwater cutting,  
Power control to maintain a constant power level during coagulation and devitalization,  
It has a cabinet for a bottle of argon, a rail for mounting additional equipment and 2 baskets for accessories,  
Diathermy has the following operating modes:  
AUTO CUT,  
BIPOLAR CUT,  
ARGON AUTO CUT,  
SOFT COAG,  
SWIFT COAG,  
FORCED COAG,  
BIPOLAR SOFT COAG,  
BIPOLAR FORCED COAG,  
FORCED APC,  
PRECISE APC,  
PULSED APC,  
ARGON SOFT COAG,  
ARGON SWIFT COAG,  
ARGON FORCED COAG,  
Specification:  
Maximum cutting power: 200 W at 500 Ω,  
Maximum coagulation power: 120 W at 500 Ω,  
Maximum bipolar power: 120 W at 75 Ω,  
HF frequency: 350 kHz,  
Safety system: NESSY,  
Mobile, on castors, 2 of which with brake,  
Included:  
Factory new monopolar electrode cable with ERGO PEN BOWA holder - two buttons, 4 mm shaft (no electrodes included, can be retrofitted at additional cost),  
Brand-new cable with BOWA clip for disposable neutral electrodes (no electrodes included, can be retrofitted at additional cost),  
Brand-new bipolar cable, BOWA forceps (bipolar tweezers can be retrofitted at additional cost),  
ERBE argon gas regulator with pressure gauge and gas hose with quick coupling,  
Argon cylinder,  
Foot attachment,  
Instruction manual in English (PDF),  
Power cord,  
Dimensions: 63 x 60 x 118 cm,  
Weight: 60.5 kg,  
Comes with a current Technical Passport issued,  
Warranty:  
6 months for the domestic market - Poland,  
1 month for the international market,  
Financing options (Poland only): Installments, Leasing, Loan,  
If you have any questions, please do not hesitate to contact us!  
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.  
WE SPEAK ENGLISH  
HABLAMOS ESPANOL  
WIR SPRECHEN DEUTSCH  
ON PARLE FRANCAIS  
WIJ SPREKEN NEDERLANDS  
МЫ ГОВОРИМ ПО РУССКИ  
МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ  
WE SHIP WORLDWIDE

