



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67
E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006
kw@praiston.pl



Infusion pump FRESENIUS PILOTE C IEC (Reconditioned)

Used,
Technical condition: very good,
Visual condition: good,
French production,
Power supply: 230 V,
Emergency power supply: through 1 new built-in 6 V / 1.3 Ah battery,
Frequency: 50/60 Hz,
Power: 23 VA,
The Pilote C IEC is a syringe infusion pump for precise dosing of fluids,
Syringes supported:
20 mL: BD PLASTIPAK, BRAUN OMNIFIX, BRAUN PERFUSOR, MONOJECT, TERUMO,
50 mL: BD PLASTIPAK, BD PERFUSION, BRAUN OMNIFIX, BRAUN PERFUSOR, MONOJECT, TERUMO,
Acoustic and visual alarms, e.g. occlusion pressure alarm, volume limit, KVO,
Equipped with syringe size reader,
Anti-bolus system,
Comfortable carrying handle and RS-232 port,
Easy to read digital display,
French keyboard on control panel,
Control panel with visible signs of use in the form of damaged foiling,
Specifications:
Infusion rate: 0.1 - 200 mL / h,
Bolus: 0.1 - 800 mL / h,
Infusion volume: 0.1 - 999 mL,
KVO: 1 mL / h,
Infusion rate accuracy: $\pm 1\%$ on drive, $\pm 2\%$ on syringe,
Pause: 1 min - 9 h 59 min,
Dimensions: 32.5 x 19 x 12 cm,
Weight: 3.5 kg,
Cat. 03,
Has a current Technical Passport issued,
Warranty:
6 months for the domestic market - Poland,
1 month for the international market,
Financing options (Poland only): Installments, Leasing, Loan,
If you have any questions, please do not hesitate to contact us!
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.
WE SPEAK ENGLISH
HABLAMOS ESPANOL
WIR SPRECHEN DEUTSCH
ON PARLE FRANCAIS
WIJ SPREKEN NEDERLANDS
МЫ ГОВОРИМ ПО РУССКИ
МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ
WE SHIP WORLDWIDE

