



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



PRAISTECH FIBI medical cabinet (New)

New,
Made in Poland,
Fibi medicine cabinet has durable, metal construction,
Equipped with double glazed doors with a lock,
Four shelves provide ample storage space,
Shelves can be adjusted in height thanks to perforated rails with holes every 1,5 cm,
2 keys included,

Dimensions:

Shelves: 90 x 37.5 cm,

General: 90 x 40 x 185 cm,

Weight: 44.5 kg,

Holds a current Technical Passport,

Warranty: 24 months,

Finance options (Poland only): Instalments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

