



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67  
E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006  
[kw@praiston.pl](mailto:kw@praiston.pl)



## **Video processor and light source OLYMPUS OTV-SI (Reconditioned)**

Used,  
Technical condition: very good,  
Visual condition: very good,  
Japanese production,  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Olympus OTV-SI is a device that integrates an image processor and a light source to save cabinet space,  
Intuitive operation,  
Main and emergency lighting: 2 halogen bulbs 15 V / 150 W,  
Suitable for flexible and rigid endoscopes,  
Has the following functions:  
Operating modes: manual and automatic,  
Contour enhancement for easier observation of tissue,  
White balance,  
Gradual adjustment of light intensity,  
It has the following inputs and outputs:  
2 x Y/C,  
VIDEO OUT,  
SYSTEM,  
H. COPY,  
VTR REMOTE,  
Socket for keyboard connection,  
Keypad on control panel in language: English,  
Included:  
OLYMPUS Model MAJ-1124 keyboard with silicone cover,  
Power cord,  
Dimensions: 43 x 30 x 17 cm,  
Set weight: 13 kg,  
It has a current Technical Passport issued,  
Warranty:  
6 months for the domestic market - Poland,  
1 month for the international market,  
Financing options (Poland only): instalments, leasing, loan,  
If you have any questions, please do not hesitate to contact us!  
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.  
WE SPEAK ENGLISH  
HABLAMOS ESPANOL  
WIR SPRECHEN DEUTSCH  
ON PARLE FRANCAIS  
WIJ SPREKEN NEDERLANDS  
МЫ ГОВОРИМ ПО РУССКИ  
МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ  
WE SHIP WORLDWIDE

