



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Certified ZHENRUI COVID-19 antigen test(s) (cassette) (1 pc) (New)

New,
P.R.C. production,
Zhenrui tests are for the rapid qualitative detection of Covid-19 antigen in a nasal, nasopharyngeal or saliva sample,
Designed for use in professional medical facilities,
The result can be read after 20 minutes,
Specification:
Sensitivity: 96 %,
Specificity: 97 %,
Accuracy: 96.5 %,
Included:
Test cassette,
Nasal swab for taking a sample from the nasopharyngeal cavity,
Swab stick for taking a sample from the nasopharyngeal cavity,
Collector to take saliva sample,
Squeezer,
Extraction tube with buffer,
Instruction manual in languages: PL, EN, DE, FR, IT, ES, PT, NL
If you have any questions, please do not hesitate to contact us!
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.
WE SPEAK ENGLISH
HABLAMOS ESPANOL
WIR SPRECHEN DEUTSCH
ON PARLE FRANCAIS
WIJ SPREKEN NEDERLANDS
МЫ ГОВОРИМ ПО РУССКИ
МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ
WE SHIP WORLDWIDE

