



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **OLYMPUS sterilisation container (Reconditioned)**

Used,  
Technical condition: very good,  
Visual condition: very good,  
Made in Japan,  
Olympus container is used for sterilisation and storage of rigid endoscopes,  
Equipped with rubber holders for mounting the endoscope,  
4 clips secure the container against accidental opening,  
The container has signs of use,  
Dimensions: 67 x 21 x 7 cm,  
Weight: 1 kg,  
Comes with a current Technical Passport issued,  
Warranty:

6 months for the domestic market - Poland,  
1 month for the international market,

Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPAÑOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANÇAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

