



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



ACL 19" ACL OR-PC 19LP medical computer (Reconditioned)

Used,
Technical condition: very good,
Visual condition: very good,
German production,
Power supply: 230 V,
Frequency: 50/60 Hz,
ACL is a medical, fanless, compact PC,
It is designed for continuous 24/7 operation in a patient environment,
Suitable for convection cooling,
Equipped with an internal 19" wide viewing angle TFT display,
Has an easy to clean design,
Has the following inputs and outputs:
HDMI,
DVI,
LAN - 2 pcs,
USB 2.0 - 4 pcs,
COM 1 and 2,
SOUND,
Specification:
Resolution: 1280 x 1024 px,
Contrast ratio 1500:1,
Brightness 300 cd/m²,
Viewing angle: 178°/178°,
Processor: IntelCore i7,
RAM: 4 GB,
System: Windows 10 Pro,
The computer case has minor signs of use,
Keyboard on the control panel in English,
Included power cord,
Dimensions: 43 x 36 x 10 cm,
Weight: 10 kg,
Has a current Technical Passport issued,
Warranty:
6 months for the domestic market - Poland,
1 month for the international market,
Financing options (Poland only): Installments, Leasing, Loan,
If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

