



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



VALLEYLAB diathermy trolley (Reconditioned)

Used,
Technical condition: very good,
Visual condition: very good,
U.S.A. production,
Valleylab is a functional trolley for storing and moving surgical diathermy,
Equipped with:
Upper table top with a working surface of 49.5 x 44 cm,
Mounting holes,
A shelf for cables and accessories of 11 cm in height,
Lower table top 54 x 46 cm,
Handles for organising cables and moving the trolley,
The lower top has walls to protect the items placed on it from falling,
Mobile on castors, 2 of which have brakes,
Dimensions: 67 x 49.5 x 89 cm,
Weight: 19 kg,
It has a current Technical Passport issued,
Warranty:

6 months for domestic market - Poland,

1 month for the international market,

Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

