



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Treatment lamp MAQUET LUCEA 10 (ceiling) (Reconditioned)

Used,

Technical condition: very good,

Visual condition: very good,

Swedish production,

Power supply: 230 V,

Frequency: 50/60 Hz,

Power: 22 VA,

Maquet Lucea 10 is a ceiling treatment lamp,

Thanks to its good colour reproduction it is particularly dedicated to dermatology,

The use of LEDs means that the light does not emit heat, which helps to distinguish between healthy and varicose veins,

It is possible to set the bowl in any position thanks to adjustment in many planes,

It has a two-part arm with a total length of 80 cm,

The bowl is equipped with 10 LEDs producing the maximum light intensity of 50,000 lx,

The lamp canopy bears minor signs of use,

Specifications:

Colour temperature: approx. 4500 K,

Colour rendering index: 95 Ra,

Light field: Ø 13 cm,

Included:

Column for hanging the lamp,

Grille,

Assembly kit,

Dimensions:

Canopy: 22 x 17 cm,

General: 72 x 26 x 15 cm,

Weight: 3 kg,

Comes with a current Technical Passport issued,

Warranty:

6 months for domestic market - Poland,

1 month for the international market,

Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

