



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Certified ZHUHAI LITUO COVID-19 antigen test(s) (colloidal gold) (5 pcs.) (New)

New,
P.R.C. production,
Zhu Hai is a diagnostic colloidal gold cassette test with CE certification for personal use,
It is designed for direct and qualitative detection of SARS-CoV-2 nucleocapsid antigens in suspected COVID-19 patients,
The sample is obtained by swabbing the front of the nose,
Reading of the result possible after 10 - 15 minutes,
Easy and self-sampling reduces the risk of virus transmission during the test,

Specification:

Sensitivity: 94.53 %,

Specificity: 99 %,

Included:

Test cassette - 5 pcs,

Swab stick - 5 pcs,

Extraction tube with buffer, stopper and dropper - 5 pcs,

Leaflet in German,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

