



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Patient monitor cardiomonitor MED CHOICE MMED6000DP-M10 (Reconditioned)

Used,

Technical condition: very good,

Visual condition: very good,

P.R.C. production,

Power supply: 230 V or through 2 new replaceable 12 V / 2.4 Ah batteries (one included, can be retrofitted at extra cost),

Frequency: 50/60 Hz,

Power: 110 VA,

Med Choice is a portable patient monitor, which thanks to its handiness and reliability is widely used in first aid stations, clinics, hospitals, medical organisations, among others,

Equipped with the following modules with full wiring and sensors:

ECG - 3 - lead,

SPO2 - Factory new silicon sensor,

NIBP - Cuff with air hose,

TEMP 1/2 - No sensor included (can be retrofitted at additional cost),

Audible and visual alarms,

FREEZE function,

Has 2 USB 2.0 ports and AUX and VGA sockets,

Wired network connection possible,

It has a 10.6" TFT screen,

English language menu,

Keypad on the control panel in: English,

Overall dimensions: 31.5 x 18 x 30 cm,

Set weight: 5 kg,

Cat. 01,

It has a current Technical Passport issued,

Warranty:

6 months for the domestic market - Poland,

1 month for the international market,

Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

