



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67
E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006
kw@praiston.pl



Video processor and light source OLYMPUS OTV-SI (Reconditioned)

Used,
Technical condition: very good,
Visual condition: very good,
Japanese production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Olympus OTV-SI is a device that integrates an image processor and a light source to save cabinet space,
Intuitive operation,
Main and emergency lighting: 2 halogen bulbs 15 V / 150 W,
Suitable for flexible and rigid endoscopes,
Has the following functions:
Operating modes: manual and automatic,
Contour enhancement for easier observation of tissue,
White balance,
Gradual adjustment of light intensity,
It has the following inputs and outputs:
2 x Y/C,
VIDEO OUT,
SYSTEM,
H. COPY,
VTR REMOTE,
Socket for keyboard connection,
Keyboard on the control panel in English,
The unit has minor signs of use in the form of scratches on the cover,
Included:
OLYMPUS Keyboard Model MAJ-1124,
Power cord,
Dimensions: 43 x 30 x 17 cm,
Set weight: 13 kg,
Has a current Technical Passport issued,
Warranty:
6 months for the domestic market - Poland,
1 month for the international market,
Financing options (Poland only): instalments, leasing, loan,
If you have any questions, please do not hesitate to contact us!
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.
WE SPEAK ENGLISH
HABLAMOS ESPANOL
WIR SPRECHEN DEUTSCH
ON PARLE FRANCAIS
WIJ SPREKEN NEDERLANDS
МЫ ГОВОРИМ ПО РУССКИ
МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ
WE SHIP WORLDWIDE

