



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



AIRCAST high foot and knee brace (Reconditioned)

Used,

Technical condition: very good,

Visual condition: very good,

Made in Mexico,

The Aircast foot and knee brace is an excellent alternative to plaster immobilisation,

Four-chamber air system providing massaging compression to reduce swelling faster,

Fits right and left leg in most adult patients,

Foam lining and air chambers for comfortable walking and effective swelling reduction,

Built-in pump to adjust air chamber compression,

Durable, semi-rigid shoe construction for optimal stability

SoftStrike technology (air-filled sole) cushions and dissipates shock while walking,

Rubber, lightweight cradle sole prevents slipping and forces a natural gait,

Size S: 35.5 - 39.5 cm,

Dimensions: 27 x 13 x 40 cm,

Weight: 1 kg,

Comes with a current Technical Passport issued,

Warranty:

6 months for domestic market - Poland,

1 month for the international market,

Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

