



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Headrest ALVO 4-09-040 (Reconditioned)

Used,

Technical condition: very good,

Visual condition: very good,

Made in Poland,

Alvo is a specialist headrest for neurological, laryngological and ophthalmological surgeries,

Anatomically shaped elements,

Manual adjustment of the headrest position in a wide range:

Angle of inclination,

Height,

Distance from the table,

Adjustable spacing of headrest elements for a better fit to the patient's head,

Included adapter for securing the headrest,

Dimensions: 50 x 37 x 17 cm,

Weight: 6.5 kg,

Comes with a current Technical Passport issued,

Warranty:

6 months for the domestic market - Poland,

1 month for the international market,

Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

