



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67
E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006
kw@praiston.pl



Cardiac patient monitor WELCH ALLYN PROPAQ CS 242 (Reconditioned)

Used,
Technical condition: very good,
Visual condition: very good,
U.S.A. production,
Power supply: 230 V,
Emergency power supply: 2 new built-in 8 V / 2.7 Ah batteries,
Frequency: 50/60 Hz,
The Propaq CS 242 is a compact patient monitor that allows you to control many of the patient's vital signs even under extreme conditions,
The device has:
8.5" LCD touchscreen display with two small scratches that do not affect the comfort of the monitor,
English language menu,
Possibility of connecting to a defibrillator through the DEFIB SYNCHRO socket,
Acoustic and visual alarms,
Handle for moving and hanging,
Equipped with the following modules with full wiring and sensors:
ECG / ECG RESP - New 3 - lead cable with snap type connectors,
SPO2 - New clip-on sensor,
NIBP - New cuff with air hose,
TEMP 1/2 - No sensor included (can be retrofitted at additional cost),
The button that activates the device has a damaged protective rubber band,
Included:
English user manual (PDF),
Power supply 16 - 24 V,
Power cord,
Dimensions: 24 x 20 x 20 cm,
Weight: 6 kg,
Comes with a current Technical Passport issued,
Warranty:
6 months for the domestic market - Poland,
1 month for the international market,
Financing options (Poland only): Installments, Leasing, Loan,
If you have any questions, please do not hesitate to contact us!
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.
WE SPEAK ENGLISH
HABLAMOS ESPANOL
WIR SPRECHEN DEUTSCH
ON PARLE FRANCAIS
WIJ SPREKEN NEDERLANDS
МЫ ГОВОРИМ ПО РУССКИ
МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ
WE SHIP WORLDWIDE

