



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67
E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006
kw@praiston.pl



Rehabilitation bed STIEGELMEYER 5085-04-85-00-0 (Reconditioned)

Used,
Technical condition: very good,
Visual condition: good,
German production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Stiegmeyer is a comfortable rehabilitation bed with electric height adjustment,
Two-section,
Equipped with wooden panels on the head and leg side with cut-outs for comfortable bed movement,
Sturdy and stable construction ensures safety and comfort,
Folding handrails on both sides of the bed,
Bumpers in each corner,
Continuous electric adjustment with a wired remote control:
Height of the bed in the range of 40 - 80 cm,
Back section angle 68°,
The bed can be fitted with a mattress and a patient arm for an additional charge,
The bed and the remote control have visible signs of use,
Mobile, on 4 wheels with two central brakes,
Dimensions:
Bed: 193 x 90 cm,
Overall: 209 x 102 x 90/130 cm,
Weight: 126 kg,
Comes with a current Technical Passport issued,
Warranty:
6 months for domestic market - Poland,
1 month for the international market,
Financing options (Poland only): Installments, Leasing, Loan,
If you have any questions, please do not hesitate to contact us!
If you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.
WE SPEAK ENGLISH
HABLAMOS ESPANOL
WIR SPRECHEN DEUTSCH
ON PARLE FRANCAIS
WIJ SPREKEN NEDERLANDS
МЫ ГОВОРИМ ПО РУССКИ
МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ
WE SHIP WORLDWIDE

