



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Rescue stretcher FERNO S-265 (Reconditioned)

Used,
Technical condition: very good,
Visual condition: good,
English production,
Very light thanks to the frame made of aluminium,
Can be separated into two parts,
Connection of the upper part of the stretcher with the lower part through hinges, allowing to fold the stretcher in half for transportation,
Four-stage length adjustment within the range 165 - 200 cm,
The stretcher has visible signs of use,
Includes instructions for use in English (PDF),
Dimensions:
General: 164 x 42 x 6 cm,
Folded length: 120 cm,
Weight: 8.5 kg,
Comes with a current Technical Passport issued,
Warranty:
6 months for domestic market - Poland,
1 month for the international market,
Financing options (Poland only): Installments, Leasing, Loan,
If you have any questions, please do not hesitate to contact us!
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.
WE SPEAK ENGLISH
HABLAMOS ESPANOL
WIR SPRECHEN DEUTSCH
ON PARLE FRANCAIS
WIJ SPREKEN NEDERLANDS
МЫ ГОВОРИМ ПО РУССКИ
МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ
WE SHIP WORLDWIDE

