



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



**TOSHIBA SSA-790A APLIO XG ultrasound system (Reconditioned)**

Used,  
Technical condition: very good,  
Visual condition: good,  
Japanese production, (2011 r.),  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Power: 1500 VA,  
The Aplio XG SSA-790A provides excellent clinical results supported by a wide range of unique, clinically validated technologies,  
A full range of imaging features enables fine details of tissue and vascular structures to be precisely visualised for faster and more accurate diagnosis.  
Equipped with:  
19.2" LCD monitor,  
Additional 10.4" touch screen display,  
5 active head slots,  
Retractable alphanumeric keyboard,  
DVD drive,  
2 USB 2.0 ports  
Built-in speakers,  
Drawer with dimensions of 30 x 9.5 x 10 cm,  
Transport handle,  
Features:  
2D,  
M,  
Pulsed Doppler (PW),  
Continuous wave (CW),  
Colour flow imaging (CDI),  
Power,  
Tissue Doppler (TDI),  
Harmonic Imaging (THI),  
Automatic image optimisation (QuickScan),  
Wired network connectivity,  
Menu language: English,  
Keypad on control panel in English,  
Visible signs of use,  
Mobile, on 4 wheels with central brake and directional drive capability,  
Included:  
MITSUBISHI Videoprinter Model P95,  
Phase head TOSHIBA Model PST-30BT:  
Frequency range:2 - 4.8 MHz,  
Transcranial and cardiac imaging,  
TOSHIBA Model PLT-704SBT Linear Head:  
Frequency range:4.8 - 11.0 MHz,  
Vascular imaging,  
TOSHIBA Model PVT-661VT endovaginal transducer:  
Frequency range:3.6 - 8.8 MHz,  
Endovaginal and endorectal imaging,  
TOSHIBA Model PLT-1204BT (18L7) linear transducer with compromised coating:  
Frequency range:7.2 - 14 MHz,  
Vascular and small area imaging,  
The TOSHIBA Model PLT-805AT (12L5) transducer shown on the images is not included,  
Dimensions: 112 x 56 x 143/167 cm,  
Weight of set: 193 kg,  
Has a current Technical Passport issued,  
Warranty:  
6 months for the domestic market - Poland,  
1 month for the international market,  
Financing options (Poland only): instalments, leasing, loan,  
If you have any questions, please do not hesitate to contact us!  
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.  
WE SPEAK ENGLISH  
HABLAMOS ESPANOL  
WIR SPRECHEN DEUTSCH  
ON PARLE FRANCAIS  
WIJ SPREKEN NEDERLANDS  
МЫ ГОВОРИМ ПО РУССКИ  
МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ  
WE SHIP WORLDWIDE

