



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Neuromuscular conductivity monitor TOF-WATCH (Reconditioned)

Used,
Technical condition: very good,
Visual condition: very good,
Dutch production,
Power supply: through a new built-in 9 V battery,
Frequency: 50/60 Hz,
TOF-Watch is a microprocessor-controlled, easy-to-use, portable neuromuscular conductance monitor using acceleromyography principles,
It provides precise information about the degree of relaxation,
Enables faster recovery of muscle and respiratory function,
Prevents the development of persistent neuromuscular block,
Automatically switches off after 2 hours of inactivity,
It has an easy to read 3" display with a scratch that does not affect the comfort of use,
Stimulation modes:
Quadruple series (TOF),
Power Tach Number (PTC),
Single Pulse (ST) 1 or 0.1 Hz,
Double pulse (DBS) 3.3 or 3.2,
Tetanic stimulation (TET) 50 or 100 Hz,
Specifications:
Stimulation current: 0 - 60 mA at an impedance < 5 kOhm,
Single phase pulse duration: 200 µs,
Included:
English user manual (PDF),
Patient cable,
Dimensions: 20 x 9.5 x 5 cm,
Weight: 1 kg,
Cat. 02,
Has a current Technical Passport issued,
Warranty:
6 months for the domestic market - Poland,
1 month for the international market,
Financing options (Poland only): Installments, Leasing, Loan,
If you have any questions, please do not hesitate to contact us!
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.
WE SPEAK ENGLISH
HABLAMOS ESPANOL
WIR SPRECHEN DEUTSCH
ON PARLE FRANCAIS
WIJ SPREKEN NEDERLANDS
МЫ ГОВОРИМ ПО РУССКИ
МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ
WE SHIP WORLDWIDE

