



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



LED light source with remote control HCM MEDICA HCM9698 (Reconditioned)

Used,
Technical condition: very good,
Visual condition: very good,
Produced by P.R.C.,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 120 W,
HCM Medica is a light source with remote control,
It is an excellent alternative to light sources based on xenon or halogen lighting,
Possibility to change the light intensity without approaching the device,
High brightness comparable to a 350 W xenon light source,
Does not emit UV or IR radiation,
Keypad on the control panel in language: English,

Specifications:

Lifetime: 60000 h,
Colour temperature: 5000 - 6000 K,
Intensity: 110,400 lux,
Luminous flux: 2484 Lm,

Included:

Wireless remote control,
Instruction manual in Polish,
Power supply cable,
Dimensions: 36 x 36 x 14 cm,
Weight: 6.5 kg,
Comes with a current Technical Passport issued,
Warranty:

6 months for domestic market - Poland,

1 month for the international market,

Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE



