



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Wrist support STRYKER 1068-250 (Reconditioned)

Used,

Technical condition: very good,

Visual condition: very good,

Production U.S.A,

Wrist support is used as a support for surgeon's forearms working in standing or sitting position behind patient's head,

Used during ophthalmic and head and neck surgical procedures,

Half-ring shape provides a wide range of support,

Infinitely variable manual adjustment of the distance between the support and the patient

Compatible with the head section of most Stryker eye, head and neck stretchers,

The support leg has visible signs of use,

Dimensions: 41 x 34 x 19 cm,

Weight: 2 kg,

Has a current Technical Passport issued,

Warranty:

6 months for domestic market - Poland,

1 month for the international market,

Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

