



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **LEISEGANG 1H colposcope (Reconditioned)**

Used,  
Technical condition: very good,  
Visual condition: very good,  
German production,  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Power: 40 VA,  
The Leisegang 1H is a chair-mounted colposcope for gynaecological examinations,  
It has a green filter,  
Mounted on a 3-part pantograph arm with a total length of approx. 80 cm,  
Light source: 12 V / 20 W halogen bulb,  
Equipped with a handle for better control of the device,  
Continuous manual adjustment:

The height of the colposcope head on the stick by means of a knob and a handle within a range of 33 cm,

Light intensity,

Pupil distance and eyepiece correction,

One of the eyepieces has a crack,

Dimensions: 58 x 23 x 73/106 cm,

Weight: 13 kg,

Cat. 04,

Has a current Technical Passport issued,

Warranty:

6 months for the domestic market - Poland,

1 month for the international market,

Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

