



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MEDELA SYMPHONY breast pump (Reconditioned)

Used,
Technical condition: very good,
Visual condition: very good,
Swiss production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Emergency power supply: through 2 built-in rechargeable batteries 6 V 1.2 Ah (not included, can be retrofitted at additional cost),
Medela Symphony was created to support mothers during breastfeeding: it helps to initiate, develop and maintain lactation at the right level,
Perfect for prolonged and frequent pumping,
Uses initiation technology and 2-Phase Expression to mirror the baby's natural sucking rhythm and help stimulate milk flow,
Overflow protection,
Allows double and single pumping,
Menu in Polish,
Comes with a cord holder,
Negative pressure generated: (-50) - (-250) mmHg,
Mobile, on a mobile base on 4 double castors, 2 of which with brake,
Included:
Symphony PLUS programme card,
Bottle holder,
Operating instructions in several languages (PDF),
Power cord,
Dimensions:
Breast pump: 21 x 14 x 25 cm,
General: 37 x 37 x 92 cm,
Weight: 6 kg,
Comes with a current Technical Passport issued,
Warranty:
6 months for domestic market - Poland,
1 month for the international market,
Financing options (Poland only): Installments, Leasing, Loan,
If you have any questions, please do not hesitate to contact us!
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.
WE SPEAK ENGLISH
HABLAMOS ESPANOL
WIR SPRECHEN DEUTSCH
ON PARLE FRANCAIS
WIJ SPREKEN NEDERLANDS
МЫ ГОВОРИМ ПО РУССКИ
МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ
WE SHIP WORLDWIDE

