



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Patient hand strap SCHMITZ 101.121 (Reconditioned)

Used,
Technical condition: very good,
Visual condition: very good,
Made in Germany,
The patient hand belt is used for immobilisation of the upper limb at wrist level,
It consists of a belt equipped with Velcro and an assembly holder with a fastener,
The holder can be removed and the belt washed at 95° C,
Infinitely adjustable clamp,
Belt made of polyester and polyamide,
Dimensions:
Belt width: 11,5 cm,
Overall: 24.5 x 11.5 x 4 cm,
Weight: 0.5 kg,
Has a current Technical Passport issued,
Warranty:
6 months for domestic market - Poland,
1 month for the international market,
Financing options (Poland only): Installments, Leasing, Loan,
If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

