



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Side support MAQUET 1002.11C0 (Reconditioned)

Used,
Technical condition: very good,
Visual condition: very good,
German production,
The side support provides a safe and comfortable position for the hips and other parts of the patient's body in procedures requiring optimal lateral stability,
It has an integrated mattress with dimensions of 21.5 x 10 x 5 cm,
Gentle contouring for better pressure distribution,
Support can be rotated,
Removable pad makes it easy to keep the accessory clean,
Fitted with an 11 cm long pole (fixing not supplied),
The stand post bears slight signs of use,
Dimensions: 21.5 x 10 x 8 cm,
Weight: 1 kg,
Cat. 01,
Has a current Technical Passport issued,
Warranty:
6 months for the domestic market - Poland,
1 month for the international market,
Financing options (Poland only): Installments, Leasing, Loan,
If you have any questions, please do not hesitate to contact us!
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.
WE SPEAK ENGLISH
HABLAMOS ESPANOL
WIR SPRECHEN DEUTSCH
ON PARLE FRANCAIS
WIJ SPREKEN NEDERLANDS
МЫ ГОВОРИМ ПО РУССКИ
МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ
WE SHIP WORLDWIDE

