



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Blood collection station U.B.M ST-02 (Reconditioned)

Used,
Technical condition: very good,
Visual condition: very good,
Made in Poland,
U.B.M. stand enables efficient and easy procedure of blood sampling on both right and left forearm,
New upholstery with medical certificate in beige colour (possibility of changing the upholstery for a new one from our colour chart at extra charge),
Equipped with 2 armrests adjustable in height from 57 - 77 cm,
Armrests placed on a pillar allowing rotation in a vertical plane,
All elements are mounted on a sturdy base of powder-coated steel profiles,
Rubber feet prevent the stand from slipping,
Dimensions:
General: 80 x 74 x 81 cm,
Seats: 39 x 40 x 4.5 cm,
Armrest: 44 x 17 x 3 cm,
Maximum load capacity: 120 kg,
Weight: 17 kg,
Cat. 03,
It has a current Technical Passport issued,
Warranty:
6 months for the domestic market - Poland,
1 month for the international market,
Financing options (Poland only): Installments, Leasing, Loan,
If you have any questions, please do not hesitate to contact us!
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.
WE SPEAK ENGLISH
HABLAMOS ESPANOL
WIR SPRECHEN DEUTSCH
ON PARLE FRANCAIS
WIJ SPREKEN NEDERLANDS
МЫ ГОВОРИМ ПО РУССКИ
МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ
WE SHIP WORLDWIDE

