



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Recording system central unit MED X CHANGE DRS01 with ELO 15" monitor (Reconditioned)

Used,

Technical condition: very good,

Visual condition: very good,

U.S.A. production,

Power supply: 230 V,

Frequency: 50/60 Hz,

It has the following inputs and outputs:

Ports to connect mouse, keyboard and switches,

CMIS,

VGA,

LAN,

LINE IN,

AUDIO OUT,

MIC,

DVI OUT,

USB-B,

HDMI - 2 pcs,

USB 2.0 - 8 pieces, including 4 on the front panel,

Equipped with a 15" ELO Model ENTUITIVE LCD monitor:

Touch screen display with a single point damage that does not affect the comfort of use of the monitor,

Navigation buttons on the side of the screen,

It has the following inputs and outputs:

DVI,

VGA,

USB-B,

AUDIO,

Specifications:

Response time: 8 ms,

Processor: IntelCore 2 Duo,

RAM: 4 GB,

System: Windows XP Pro,

Menu selection in EN, DE, IT, ES, PT, SL,

The control unit has signs of use in the form of scratches on the housing,

Included:

ELO Monitor Model ENTUITIVE,

Keyboard with touchpad MED X CHANGE,

Connection cables,

12 V power supply,

Dimensions:

Monitor: 40 x 26.5 x 21 cm,

Central unit: 36 x 20 x 33.5 cm,

Weight of set: 18 kg,

Comes with a current Technical Passport issued,

Warranty:

6 months for the domestic market - Poland,

1 month for the international market,

Financing options (Poland only): instalments, leasing, loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

