



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Infusion pump (syringe) B. BRAUN PERFUSOR SECURA FT 12/230 V (Reconditioned)

Used,
Technical condition: good,
Visual condition: good,
German production,
Power supply: 12 / 230 V or through 1 built-in rechargeable battery 7.2 V / 1.2 Ah,
Frequency: 50/60 Hz,
Perfusor Secura FT is a syringe-operated, single-channel infusion pump for precise fluid dosing,
It supports B. BRAUN with a capacity of 50 ml,
It has:
Easy to read display,
Foiled control panel,
Comfortable carrying handle,
Built-in alarms and prealarms, e.g. 3 minutes before end of infusion,
Infusion rate: 0.1 - 99.9 mL / h,
Quick reference guide on the control panel in English,
The pump has visible signs of use,
Included:
User manual in languages: EN, DE (PDF),
Power cord,
Dimensions: 31.5 x 15.5 x 17 cm,
Weight: 2.5 kg,
Cat 04,
Has a current Technical Passport issued,
Warranty:
6 months for the domestic market - Poland,
1 month for the international market,
Financing options (Poland only): Installments, Leasing, Loan,
If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

