



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Carboxytherapy device FENIX (Reconditioned)

Used,
Technical condition: very good,
Visual condition: very good,
Power supply: 230 V,
Frequency: 50/60 Hz,
Fenix is a highly technologically advanced medical device for carboxytherapy, which has ready-made treatment programmes and gas heating guarantees the absence of pain sensations. Thanks to these features the treatment becomes simple and safe for both the patient and the person carrying out the procedure,
Carboxytherapy is a treatment which has already gained an established position in aesthetic medicine surgeries. It works by injecting controlled doses of pure carbon dioxide (CO2) subcutaneously or intradermally. CO2 is produced by our cells as a result of burning oxygen and is excreted by our lungs. Because it is inert to us, it can also be safely injected into the skin,
Carboxytherapy is a minimally invasive and 100% natural treatment that does not require convalescence. Scientific studies have shown that it improves skin elasticity and reactivates microcirculation locally stimulates collagen fibres and breaks down the adipocyte membrane, causing a lipolytic and lipoclastic effect, without damaging connective tissue, vessels or surrounding nervous structures,
It has 5 working modes:
4 automatic programmes,
Manual programme,
Advantages of the system:
Heated gas - guarantee of a painless treatment,
Pneumatic gun for quick and precise application of CO2,
Ready-made treatment programmes - precision and safety every time,
Control of flow and volume of carbon dioxide,
Innovation - Three Way Injection - treatment option for fatty tissue of the abdomen area,
Carboxytherapy of sensitive areas such as eyes, scalp,
Always strictly specified doses of carbon dioxide,
Synergic effect of carboxytherapy and mesotherapy in one treatment,
Injection with a pistol:
High precision of gas administration,
Easy and quick injection of carbon dioxide,
Comfort and safety for the patient and the doctor,
Micrometric regulation - full control of injection depth,
Injection of strictly specified doses of carbon dioxide in both automatic and manual mode,
Built-in manometer - control of ejaculation pressure up to 0.5 MPa,
Programs allowing from 66 to 420 injections per minute in automatic mode,
Lightweight and easy to use,
Continuous power supply,
Menu in language: English,
Specifications:
Volume range: 10 - 3000 cc,
Speed range: 5 - 100 cc/m,
Mobile, on a trolley having:
Space for gas cylinder,
4 castors, 2 of which with brake,
Included:
SHERMAN pressure gauge,
Air filter,
Foot attachment,
Pistol with integrated pressure gauge,
Connection cable with hose,
Treatment protocols,
Power cord,
Dimensions:
Apparatus: 32 x 21 x 25.5 cm,
General: 38 x 31 x 109 cm,
Weight: 51 kg,
Has a current Technical Passport issued,
Warranty:
6 months for domestic market - Poland,
1 month for the international market,
Financing options (Poland only): Installments, Leasing, Loan,
If you have any questions, please do not hesitate to contact us!
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.
WE SPEAK ENGLISH
HABLAMOS ESPANOL
WIR SPRECHEN DEUTSCH
ON PARLE FRANCAIS
WIJ SPREKEN NEDERLANDS
МЫ ГОВОРИМ ПО РУССКИ
МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ
WE SHIP WORLDWIDE

