



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Ultrasound therapy device SIEMENS SONOSTAT 733 (Reconditioned)

Used,
Technical condition: very good,
Visual condition: very good,
German production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Siemens Sonostat 733 is a compact ultrasound therapy device,
Its action can be successfully used in the treatment of inflammatory, rheumatic diseases of the musculoskeletal system, injuries, inflammatory diseases of peripheral nerves, peripheral circulatory diseases,
Intuitive operation,
Continuous or pulse mode of operation,
Continuous, manual adjustment of settings:
Power in the range 0,1 - 3 W/cm²,
Treatment time in the range of 0 - 25 min,
The set includes an ultrasound head with minor signs of use,
Dimensions: 40 x 28 x 10 cm,
Weight: 6 kg,
Comes with a current Technical Passport issued,
Warranty:
6 months for domestic market - Poland,
1 month for the international market,
Financing options (Poland only): Installments, Leasing, Loan,
If you have any questions, please do not hesitate to contact us!
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.
WE SPEAK ENGLISH
HABLAMOS ESPANOL
WIR SPRECHEN DEUTSCH
ON PARLE FRANCAIS
WIJ SPREKEN NEDERLANDS
МЫ ГОВОРИМ ПО РУССКИ
МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ
WE SHIP WORLDWIDE

